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Produktbeschreibung

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EN Instructions for use

Product description

Zetuvit Plus Silicone Border is a sterile self-adhesive absorbent dressing with a silicone interface. The dressing comprises a perforated silicone wound contact layer, a semipermeable polyurethane backing film and a green absorbent pad in between two soft layers. Acrylic adhesive is used to bind these layers together. The backing film is permeable to air but waterproof, which allows the patient to shower. The silicone layer facilitates the application of the dressing and promotes non-adherence to the wound. Above the silicone layer, the product features a white, hydrophilic non-woven fabric (viscose and polypropylene). The absorbent pad comprises cellulose and superabsorbent polyacrylate (SAP) for the absorption and retention of exudate and is wrapped in a hydrophilic cellulose tissue to distribute the liquid. It is covered with a hydrophobic polypropylene non-woven fabric to indicate the site facing away from the wound.

Properties and mode of action
Zetuvit Plus Silicone Border absorbs and retains up to high levels of exudate in the absorbent pad and the dressing is able to redistribute pressure. It features a micro-adhesive silicone interface and adhesive borders, so no additional materials are required to secure the dressing. The backing film is bacteria and showerproof. The dressing changes can be carried out atraumatically and almost painlessly.

Intended purpose
Zetuvit Plus Silicone Border is a sterile single-use self-adhesive absorbent dressing with a silicone interface for long-term treatment of acute and chronic exuding wounds. It is used on adults in hospital and at home. The dressing is suitable for use in a wide range of contexts of compression therapy. Zetuvit Plus Silicone Border can be used by healthcare professionals and lay persons under the supervision of healthcare professionals.

Indications
Zetuvit Plus Silicone Border is suitable for the treatment of a wide range of exuding wounds: acute wounds (burns, post-op, trauma wounds) and chronic wounds (venous ulcers, pressure ulcers, diabetic foot ulcers, tumour wounds).

Mode of application

- Cleanse the wound according to clinical practice. Make sure that the surrounding skin is dry and clean to ensure that the dressing sticks.
- Select the dressing's size and shape, so that the green absorbent pad exceeds the wound edges by approx. 1 to 2 cm.
- Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Remove the protective sheets according to the application instructions. Apply the adhesive side to the wound with the green side facing away.
- Press down the adhesive edges firmly and avoid wrinkles in the border to ensure the dressing stays securely in place.
- No additional fixation is required. If indicated, compression bandages can be applied as instructed by a clinician.
- The dressing can remain on the wound for up to 7 days, depending on the condition of the wound and the clinician's instructions.
- The dressing should not be used if clinically indicated or when exudate approached the rim of the green absorbent pad.
- Eine zusätzliche Fixierung ist nicht erforderlich. Bei Bedarf können gemäß ärztlicher Anweisung Kompressionsbinden angelegt werden.

Contraindications

Do not use Zetuvit Plus Silicone Border in case of hypersensitivity to any of its components.

Special precautions

- Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Before treating a wound that is not healing as expected, consult a clinician. Treatment with Zetuvit Plus Silicone Border cannot replace a causal treatment of the wound-healing impairment. If there are clinical signs of infection, the infection needs to be controlled with appropriate treatment before the dressing can be used.
- In all cases, follow the clinician's instructions.
- In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used with caution and following a clinician's recommendation.
- It is not recommended to stand on the dressing.
- Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.
- Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Side effects
If any of the following side effects occur, consult a clinician: increased exudate, maceration, increased pain, critical colonization, adherence of the dressing to the wound, itching, burning, oedema, multiple ulceration, skin tears, blisters, erythema.

- Incident reporting**
If you have any information regarding a suspected adverse event, please report it to the manufacturer and/or its authorized representative and to your national authority.
- Product disposal**
To minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable local and national laws, regulations and infection prevention standards. Dispose of the medical device with household waste.
- Date of revision of the text: 2023-09-28

AU – PAUL HARTMANN Pty. Ltd. - Macquarie Park, NSW 2113
PAUL HARTMANN Pty. Ltd. - Heywood Distribution Park, Pilsworth Road, OL10 2TT
US – HARTMANN USA Inc. - Rock Hill, SC 29730
ZA – HARTMANN South Africa - Northridge, Johannesburg, Gauteng
Internationale Behörden.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices): if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal
To minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable local and national laws, regulations and infection prevention standards. Dispose of the medical device with household waste.

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NL – PAUL HARTMANN B.V. - 6546 BC Nijmegen
BE – N.V. PAUL HARTMANN S.A. - 1480 Saintes/Sint-Reeelde
IT – PAUL HARTMANN S.p.A. - 37139 Verona